Refer Patients with Afib for Heart Surgery?



Post–Ablation Considerations:

- Not all surgical ablation is a full Maze lesion set (bi–atrial lesion set), the patient may instead have:
 - PVI lesions
 - Box lesions
 - LA lesion set
- Each may vary in efficacy, with the full Maze being the most efficacious.

The "blanking period": The first 90 days post-op where recurrence of atrial tachyarrhythmias could occur due to temporary inflammatory and proarrhythmic changes.

- Afib that occurs after the blanking period may not resolve, and those patients need further management and/or intervention.
- There is currently a lack of evidence for or against the merits of anticoagulation following surgical ablation. HRS recommends making this decision based on the patients overall stroke risk, regardless of the presence of sinus rhythm.
- Electrophysiologists have consensus on monitoring protocols, management, and follow up for post-ablation patients, which can include those concomitant surgical ablation patients.

"The most successful programs in the future might be those that employ an interdisciplinary, collaborative team approach to the treatment of AF, resulting in higher success rates for patients. Many of these patients are well read and mobile and will seek out such centers, thus increasing both catheter and surgical volumes."²

Help Patients Live Better. Longer.

Increasingly more data show that surgical ablation during heart surgery reduces mortality, risk of stroke, and other post-surgical complications. **Concomitant Surgical Ablation (SA) has a Class I Recommendation from both The Society of Thoracic Surgeons (STS) and Heart Rhythm Society (HRS)**.¹²

- Patients who undergo SA treatment may actually have reduced hospital LOS.³
- One year after CABG with SA for Afib, survival improves by 42%.
- Ten years after CABG surgery with SA survival improves by 20%.⁵
- Concomitant SA gives patients with non-paroxysmal Afib the highest chance of restored sinus rhythm.⁵⁷
- Patients with a surgically restored NSR have improved quality of life.*

Restoring sinus rhythm with surgical ablation in patients with advanced heart failure can dramatically improve heart function by reducing the burden of Afib.⁹⁻¹

During surgical ablation, the surgeon also has the ability to fully exclude and electrically isolate the left atrial appendage as a part of the ablation procedure.

Lesion Set Options.

Cox Maze IV yields the highest efficacy for Afib treatment, but literature shows progressive efficacy for each additive lesion set of the Cox Maze IV.

Approach	Reported Experiences w/ Surgical Ablation	Ablation Duration
Pulmonary Vein Isolation (PVI) PAF ~50-90% ^{1,12,13} nPAF ~60% ^{1,14}	Note: + = Time	
	nPAF ~ 60% ^{1,14}	+
Box Set Lesion (Box)	nPAF ~55–70% ^{15,16}	++
Left Atrial Lesion Set (LAL)	nPAF ~73-86% ¹⁷⁻¹⁹	+++
Bi-Atrial Lesion Set (Maze)	nPAF ~80-90% ²⁰⁻²²	++++

LAA exclusion has always been a part of the Maze procedure.

Left Atrial Appendage Management (LAAM)	Effectiveness of LAAM Modalities	
LAAM is often part of surgical ablation procedures	Epicardial Clip Exclusion: 97% (93-100%) ²³⁻³³ Excision: 74% (45-100%) successful closure ³³⁻³⁶ Staple Ligation: 56% (0-71%) successful closure ^{34,35,37} Suture Ligation: 36% (23-49%) successful closure ^{34,35,37,38}	Reported Experiences 1–5 year retro and prospective peer- reviewed publications both on and off AADs

A meta-analysis for endocardial PVI outcomes for paroxysmal Afib show about 70% success in restoring sinus rhythm.³⁹ For persistent and long-standing persistent Afib, the success rates per a meta-analysis drop to about 50% and may require multiple ablation procedures.³⁹

Less than 10% of CABG patients with Afib get concomitant surgical ablation.^{40,41} dfib receive the Class I recommended concomitant surgical ablation.^{40,42}

Collaborate with the surgeon upon referral for heart surgery to determine if the patient has Afib. Ask if surgical ablation, including LAAM, can be included during the primary heart surgery.

References:

- ¹Badhwar, V., Rankin, J.S., Damiano, R.J. Jr, Gillinov, A.M., Bakaeen, F.G., Edgerton, J.R., Philpott, J.M., McCarthy, P.M., Bolling, S.F., Roberts, H.G., Thourani, V.H., Suri, R.M., Shemin, R.J., Firestone, S., Ad, N. (2017). The Society of Thoracic Surgeons 2017 Clinical Practice Guidelines for the Surgical Treatment of Atrial Fibrillation. Ann Thorac Surg, 103(1), 329-341, DOI: 10.1016/j.athoracsur.2016.10.076.
 ²Calkins, H., et al. (2017). 2017 HRS/EHRA/ECAS/APHRS/SOLAECE Expert Consensus Statement on Catheter and Surgical Ablation of Atrial Fibrillation. Heart Rhythm, 14(10):e275-e444, DOI: 10.1016/j.hrthm.2017.05.012.
- ³ Iribarne, A. et al. (2019). Surgical Atrial Fibrillation Ablation Improves Long-Term Survival: A Multicenter Analysis. Ann of Thorac Surg, 107(1):135-42.
- ⁴Rankin, J.S., Lerner, D.J., Braid-Forbes, M.J., Ferguson, M.A., & Badhwar, V. (2017). One-year mortality and costs associated with surgical ablation for atrial fibrillation concomitant to coronary artery bypass grafting. Eur J Cardiothorac Surg, 52(3):471-7.
- ⁵Musharbash, F.N. et al. (2018). Performance of the Cox-maze IV procedure is associated with improved long-term survival in patients with atrial fibrillation undergoing cardiac surgery. J Thorac Cardiovasc Surg, 155(1):159-70.
- ⁶ Philpott, J.M. et al. (2015). The ABLATE Trial: Safety and Efficacy of Cox Maze-IV Using a Bipolar Radiofrequency Ablation System. Ann of Thorac Surg, 100(5):1541-8.
- ⁷Schill, M.R. et al. (2017). Late results of the Cox-maze IV procedure in patients undergoing coronary artery bypass grafting. J Thorac Cardiovasc Surg, 153(5):1087-94.
- *Forlani, S. et al. (2006). Conversion to Sinus Rhythm by Ablation Improves Quality of Life in Patients Submitted to Mitral Valve Surgery. Ann of Thorac Surg, 81(3):863-7.
- ⁹ Adademir, T. et al. (2019). Surgical Ablation of Atrial Fibrillation in Patients With Tachycardia-Induced Cardiomyopathy. The Annals of Thoracic Surgery, 108(2):443-50.
- ¹⁰ Marrouche, N.F. et al. (2018). Catheter Ablation for Atrial Fibrillation with Heart Failure. The New England Journal of Medicine, 378:417-27.
- ¹¹ Prabhu, S. et al. (2017). Catheter Ablation Versus Medical Rate Control in Atrial Fibrillation and Systolic Dysfunction: The CAMERA-MRI Study. Journal of the American College of Cardiology, 70(16):1953-61.
- ¹² Gillinov, A.M. et al. (2005). Surgical ablation of atrial fibrillation with bipolar radiofrequency as the primary modality. J Thorac Cardiovasc Surg, 129(6):1322-9.
- ¹³ Robertson, J.O. et al. (2013). Surgical Techniques Used for the Treatment of Atrial Fibrillation. Circ J, 77(8):1941-51.
- ¹⁴ Gillinov, A.M. et al. (2015). Surgical ablation of atrial fibrillation during mitral-valve surgery. N Engl J Med, 372(15):1399-409.
- 15 Voeller, R.K. et al. (2008). Isolating the entire posterior left atrium improves surgical outcomes after the Cox maze procedure. J Thorac Cardiovasc Surg, 135(4):870-7.
- 16 Gillinov, A.M. et al. (2006). Surgery for permanent atrial fibrillation: impact of patient factors and lesion set. Ann Thorac Surg, 82(2):502-14.
- ¹⁷ Barnett, S.D., & Ad, N. (2006). Surgical ablation as treatment for the elimination of atrial fibrillation: a meta-analysis. J Thorac Cardiovasc Surg, 131(5):1029-35.
- 18 Cox, J.L., & Ad, N. (2000). The importance of cryoablation of the coronary sinus during the Maze procedure. Semin Thorac Cardiovasc Surg, 12(1):20-4.
- ¹⁹ Ad, N., Holmes, S.D., Lamont, D., & Shuman, D.J. (2017). Left-Sided Surgical Ablation for Patients With Atrial Fibrillation Who Are Undergoing Concomitant Cardiac Surgical Procedures. Ann Thorac Surg. 103(1):58-65.
- 20 Gaynor, S.L. et al. (2015). Surgical treatment of atrial fibrillation: predictors of late recurrence. J Thorac Cardiovasc Surg, 129(1):104-11.
- ²¹Weimar, T. et al. (2011). The Cox-maze IV procedure for lone atrial fibrillation: a single center experience in 100 consecutive patients. J Interv Card Electrophysiol. 31(1):47-54.
- ²² Schill, M.R. et al. (2017). Late results of the Cox-maze IV procedure in patients undergoing coronary artery bypass grafting. J Thorac Cardiovasc Surg, 153(5):1087-94.
- ²³ Ailawadi, G. et al. (2011). Exclusion of the left atrial appendage with a novel device: early results of a multicenter trial. J Thorac Cardiovasc Surg, 142(5):1002-9.
- ²⁴ Caliskan, E. et al. (2017). Epicardial left atrial appendage AtriClip occlusion reduces the incidence of stroke in patients with atrial fibrillation undergoing cardiac surgery. Europace 20(7):e-105-e114. DOI: 10.1093/europace/eux211.
- ²⁵ van Laar, C. et al. (2018). Thoracoscopic Left Atrial Appendage Clipping: A Multicenter Cohort Analysis. JACC Clin Electrophysiol, 4(7):893-901.
- ²⁶ Ellis, C.R. et al. (2017). Angiographic Efficacy of the AtriClip Left Atrial Appendage Exclusion Device Placed by Minimally Invasive Thoracoscopic Approach, JACC Clin Electrophysiol, 3(12):1356-65.
- ²⁷ Kurfirst, V. et al. (2017). Epicardial clip occlusion of the left atrial appendage during cardiac surgery provides optimal surgical results and long-term stability. Interact Cardiovasc Thorac Surg, 25(1):37-40.
- 28 Emmert, M.Y. et al. (2014). Safe, effective and durable epicardial left atrial appendage clip occlusion in patients with atrial fibrillation undergoing cardiac surgery: first long-term results from a
- prospective device trial. Eur J Cardiothorac Surg, 45(1):126-31.
- 29 Ad, N. et al. (2015). New Approach to Exclude the Left Atrial Appendage During Minimally Invasive Cryothermic Surgical Ablation. Innovations (Phila), 10(5):323-7.
- 30 Gerdisch, M. et al, AtriClip PRO+V Left Atrial Appendage Occlusion Study. AtriCure, Inc., Post Market Field Evaluation of the PRO+V Device, PM-US-0071A-1020-G.
- ³¹ Mokracek, A. et al. (2015). Thoracoscopic Occlusion of the Left Atrial Appendage. Innovations (Phila), 10(3):179-82.
- 32 Page, S. et al. (2019). Left Atrial Appendage Exclusion Using the AtriClip Device: A Case Series. Heart Lung Circ, 28(3):430-5.
- ³³ Beaver, T.M. et al. (2016). Thoracoscopic Ablation With Appendage Ligation Versus Medical Therapy for Stroke Prevention: A Proof-of-Concept Randomized Trial. Innovations (Phila), 11(2):99-105.
 ³⁴ Kandarian, A.S., Gillinov, A.M., Pettersson, G.B., Blackstone, E., & Klein, A.L. (2008). Success of surgical left atrial appendage closure: assessment by transesophageal echocardiography. J Am Coll Cardiol. 52(11):924-9.
- 35 Cullen, M.W. et al. (2016). Left Atrial Appendage Patency at Cardioversion After Surgical Left Atrial Appendage Intervention. Ann Thorac Surg, 101:675-81.
- ³⁸ Lee, R. et al. (2016). A randomized, prospective pilot comparison of 3 atrial appendage elimination techniques: Internal ligation, stapled excision, and surgical excision. J Thorac Cardiovasc Surg, 152(4):1075-80.
- ³⁷ Healey, J.S. et al. (2005). Left Atrial Appendage Occlusion Study (LAAOS): results of a randomized controlled pilot study of left atrial appendage occlusion during coronary bypass surgery in patients at risk for stroke. Am Heart J. 150(2):288-93.
- ³⁸ Lee, R. et al. (2016). A randomized, prospective pilot comparison of 3 atrial appendage elimination techniques: Internal ligation, stapled excision, and surgical excision. J Thorac Cardiovasc Surg, 152(4):1075-80.
- ³⁸ Kirchhof, P. et al. (2016). 2016 ESC guidelines for the management of atrial fibrillation developed in collaboration with EACTS: the Task Force for the management of atrial fibrillation of the European Society of Cardiology (ESC) developed with the special contribution of the European Heart Rhythm Association (EHRA) of the ESC Endorsed by the European Stroke Organization (ESO). Europace, 18(11):1455-90.
- 40 Braid Forbes, M.J. Health Research, 2014 CMS SAF, August 2016. NIS for volume and Medicare to look back over 3 years and obtained diagnosis.
- ⁴¹McCarthy, P.M. et al. (2019). Prevalence of Atrial Fibrillation before Cardiac Surgery and Factors Associated with Concomitant Ablation. J Thorac Cardiovasc Surg, PII: S0022-5223(19)31361-3, DOI: 10.1016/J.JTCVS.2019.06.062.
- ⁴² Badhwar, V. et al. (2017). Surgical ablation of atrial fibrillation in the United States: trends and propensity matched outcomes. Ann of Thorac Surg, 104(2):493-500.

Clinical results are not predictive and individual results may vary.

ATRICURE, INC.

7555 Innovation Way Mason, Ohio 45040 Customer Service: +1 (866) 349-2342 CustomerService@AtriCure.com www.AtriCure.com

www.ACTagainstAfib.com



